



Safety Data Sheet

Rev.01 Jan.06

Page 1 (3)

1. Product identification

Product name:	<i>Copan SRK Solution 10 mL</i>
Product number:	<i>R2360</i>
Supplier:	<i>COPAN DIAGNOSTICS Inc., 26055 Jefferson Ave., 92562 Murrieta, CA, 951-696-6957</i>

2. Composition/information on ingredients

Hazardous ingredients :	Percentage on medium:	Classification :
Calcium Chloride	0,01%	R36
Sodium thioglycollate	0,0375%	R21\22

3. Hazards identification

Skin contact:	Irritation
Eye contact:	Irritation

4. First-aid measures

Skin contact:	Immediately wash with water and soap and rinse thoroughly.
Eye contact:	Rinse opened eye for 15 minutes under running water. Then consult a doctor.
Inhalation:	In case of unconsciousness place patient on side position for transportation.
Swallowing:	Immediately call a doctor.

5. Fire-fighting measures

Extinguishing media:	CO2, extinguishing powder or water spray. Fight larger fires with water or alcohol resistant foam
Protective equipment:	Wear self-contained respiratory protective device Wear fully protective suit



Safety Data Sheet

Page 2 (3)

Rev.01 Jan.06

6. Accidental release measures

Precautions: Wear protective equipment. Keep unprotected persons away

7. Handling and storage

Handling: No special measures required

Storage: Store away from oxidizing agents
Store in cool, dry conditions in well sealed receptacles.

8. Exposure controls/personal protection

Respiratory protection: N/A
Hand protection: Latex gloves
Eye protection: N/A
Skin protection: Protective work clothing (lab coat)

9. Physical and chemical properties

Appearance/odour: Characteristic
pH: 7.00 ± 0.50
Density, kg/l: Not determined
Melting point, °C: Not determined
Boiling point, °C: Not determined
Solubility: Insoluble
Flash point, °C: Not applicable

10. Stability and reactivity

Materials and conditions to avoid: No decomposition if used according to specifications
Hazardous decomposition products: No dangerous decomposition products known



Safety Data Sheet

Page 3 (3)

Rev.01 Jan.06

11. Toxicological information

Skin contact:	Irritant to skin and mucous membranes
Eye contact:	Irritating effect
Toxicological data:	Irritant
	The products shows the following dangers according to internally approved calculation methods for preparations

12. Ecological information

Ecology:	The ecological effects have not been thoroughly investigated, but currently none have been identified.
General Notes:	Not known to be hazardous to water.

13. Disposal considerations

Reccomendation:	Must adhere to state and federal regulations Must not be disposed of with solid waste
------------------------	--

14. Transport information

ADR/RID:	IMO/IMDG:	ICAO/IATA-DGR:
No	No	No

15. Regulatory information

Product related hazard informations:	The product has been classified and marked in accordance with regulations on hazardous materials
---	---

16. Other information

To the best of our knowledge, the information contained herein is accurate. However, neither Copan Italia or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user.

All materials may present unknown hazards and should be used with caution.

Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.

· Department issuing MSDS: Safety & Environment Department

· Contact: Technical Service Representative

The information in this safety data sheet has been obtained from current and reliable sources. However, the data is provided without warranty, express or implied, regarding its correctness or accuracy. It is the users' responsibility to determine safe conditions for use of this product and to assume liability for loss, injury, damage, or expense resulting from improper use of this product. EN-2141-SDS-0198